

Analytical Method for Florfenicol (Animal and Fishery Products)

1. Analytes

Florfenicol

Metabolite FFNH₂ [(1*R*,2*S*)-1-(4-methylsulfonylphenyl)-2-amino-3-fluoro-1-propanol]
(hereinafter referred to as florfenicol amine, which includes compounds that are converted by hydrolysis.)

2. Application

Animal and fishery products

3. Instrument

Liquid chromatograph-tandem mass spectrometer (LC-MS/MS)

4. Reagents

Use the reagents listed in Section 3 of the General Rules, except the following.

Sulfonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer cartridge (150 mg): A polyethylene tube of 12–13 mm inside diameter packed with 150 mg of sulfonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer or a cartridge equivalent to the specified one in separation capability.

Reference standard of florfenicol: Contains not less than 98% of florfenicol.

Reference standard of florfenicol amine: Contains not less than 97% of florfenicol amine.

5. Procedure

1) Hydrolysis and extraction

Add 20 mL of 6 mol/L hydrochloric acid to 10.0 g of sample, seal tightly, and heat at 100°C for 3 hrs. Allow to cool, add 10 mL of water and 40 mL of *n*-hexane, shake, and centrifuge at 3,000 rpm for 5 min. Discard the *n*-hexane layer, add 2 g of diatomaceous earth to the aqueous layer, mix, and filter with suction. Wash the container and the residue on the filter paper with 10 mL of acetone and water (1:1) mixed solution twice, and filter the washing with suction. Combine the filtrates, and add water to make exactly 100 mL.

2) Clean-up

i) Porous diatomaceous earth column chromatography

Take exactly a 2.5 mL aliquot of the solution obtained in 1), add 1 mL of 7.5 mol/L sodium hydroxide solution and 1 g of sodium chloride, mix well, and inject it into a porous diatomaceous earth cartridge (for 5 mL volume). After leaving the cartridge for 30 min, inject 30 mL of ethyl acetate, concentrate the eluate at below 40°C and remove the solvent. Dissolve the residue in 1 mL of 0.1 vol% acetic acid and methanol (1:1) mixed solution.

- ii) Sulfonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer column chromatography

Inject 2 mL of methanol and 3 mL of 0.1 vol% acetic acid and methanol (1:1) mixed solution into a sulfonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer cartridge (150 mg) sequentially and discard each effluent. Transfer the solution obtained in i) to the cartridge, add 3 mL of 0.1 vol% acetic acid and methanol (1:1) mixed solution and 2 mL of methanol sequentially, and discard each effluent. Then, add 5 mL of ammonia solution and methanol (1:99) mixed solution, concentrate the eluate at below 40°C, and remove the solvent. Dissolve the residue in a mixed solution of 0.1 vol% acetic acid and methanol (9:1) to make exactly 1 mL, and use this solution as the test solution.

6. Calibration curve

Prepare florfenicol amine standard solutions (0.1 vol% acetic acid and methanol [9:1]) of several concentrations, inject each solution into LC-MS/MS and make calibration curves by peak-height or peak-area method. When the test solution is prepared following the above procedure, the concentration of florfenicol amine in the test solution corresponding to 0.01 mg/kg in the sample is 0.0025 mg/L.

7. Quantification

Inject the test solution into LC-MS/MS, and calculate the concentration of florfenicol amine from the calibration curve prepared in 6.

8. Confirmation

Confirm using LC-MS/MS.

9. Measurement conditions

(Example)

Column: Octadecylsilanized silica gel, 2.1 mm inside diameter, 150 mm in length and 3 µm in particle diameter

Column temperature: 40°C

Mobile phase: A mixed solution of acetonitrile and 0.1 vol% acetic acid (1:99)

Ionization mode: ESI (+)

Major monitoring ions (*m/z*): Precursor ion 248, product ions 230, 130

Injection volume: 5 µL

Expected retention time: 5 min

10. Limit of quantification

0.01 mg/kg

11. Explanatory note

1) Outline of analytical method

The method consists of adding hydrochloric acid to the sample, heating it, hydrolyzing florfenicol and its metabolites to florfenicol amine, washing with *n*-hexane, cleaning up with a porous diatomaceous earth cartridge and a sulfonate-modified divinylbenzene-*N*-

vinylpyrrolidone copolymer cartridge and quantifying and confirming using LC-MS/MS.

2) Notes

- i) For the hydrolysis procedure, confirm that hydrolysis has been carried out sufficiently using a reference standard of florfenicol.
- ii) In hydrolysis, the solution should be heated until all sample lumps have disappeared and the color has turned dark brown to black.
- iii) If the filtration velocity is slow during suction filtration, glass fiber filter paper is recommended to use.
- iv) If bubbles are generated and the volume cannot be adjusted properly during suction filtration, centrifuge at 3,000 rpm for 5 min.
- v) Suction filtration is not required if no suspended solids are observed after washing with *n*-hexane. The procedure without suction filtration is as follows.
Discard the *n*-hexane layer, add 10 mL of acetone to the aqueous layer, and add water to make exactly 100 mL.
- vi) Mix gently by inverting when adjusting the volume since some samples may foam if shaken vigorously.
- vii) In sulphonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer column chromatography, when water remains during concentration, it is recommended to concentrate to about 0.5 mL at below 40°C, then concentrate again by adding about 2 mL of ethanol.
- viii) Sulfonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer column chromatography may be omitted when clean-up is sufficient using only porous diatomaceous earth column chromatography.
- ix) In the spike recovery test with the spiking level at MRL during the development of the analytical method, sulphonate-modified divinylbenzene-*N*-vinylpyrrolidone copolymer column chromatography was omitted.
- x) When the analytical method for florfenicol amine using LC-MS/MS was developed, the following monitoring ions were used:
for quantitative ions (*m/z*): precursor ion 248, product ion 230
for qualitative ions (*m/z*): precursor ion 248, product ion 130
- xi) The cartridge shall be washed by increasing the acetonitrile concentration in the mobile phase after eluting florfenicol amine to reduce the influence of matrix carryover in the sample in LC-MS/MS measurements.
- xii) The measurement conditions with LC-MS/MS for florfenicol used in developing the analytical method are as follows:
Column: Octadecylsilanized silica gel, 2.1 mm inside diameter, 150 mm in length and 3 µm in particle diameter
Column temperature: 40°C
Mobile phase: Initially, a mixed solution of acetonitrile and 0.1 vol% acetic acid

(1:99) for 6 min, followed by a linear gradient (1:99) to (19:1) in 6 min and hold at (19:1) for 5 min.

Ionization mode: ESI (-)

Major monitoring ions

for quantitative ions (m/z): precursor ion 356, product ion 336

for qualitative ions (m/z): precursor ion 356, product ion 185

Injection volume: 5 μ L

Expected retention time: 11 min

xiii) Food items used to develop the analytical method: cattle muscle, cattle fat, cattle liver, and eel.

12. Reference

Determination and Confirmation of Florfenicol. CLG-FLOR 1.04. Revision 04. United States Department of Agriculture, Food Safety and Inspection Service, Office of Public Health Science, 2010.

13. Type

C